



Clinical Applications

Recommended Catheters

Flushing and short term (0-1 day)

- Allows the vet or technician to flush uroliths retrograde into the bladder in the case of an obstruction
- Holes in the base plate allow suturing to the prepuce
- Luer lock on hub for easy connection to syringe or urine bag

REMEMBER TO MONITOR URINE FLOW REGULARLY

BUSTER Easy Slide Cat Catheter, open end

- ✓ Optimal rigidity for easy catheterization
- ✓ Softens and adjust to the anatomy
- ✓ Radiopaque for x-ray visibility

BUSTER Cat Catheter, open end

- ✓ Soft and flexible catheter

Medium term (1-3 days)

- Secure urine passage. Typically used in FLUTD/FIC cats, hospitalized patients and patients undergoing surgery
- The catheter has two lateral side holes and a rounded, atraumatically closed end
- Comes in various sizes with or without stylet
- Holes in the base plate allow suturing to the prepuce
- Luer lock on hub for easy connection to syringe or urine bag

BUSTER Easy Slide Cat Catheter, side holes

- ✓ Optimal rigidity for easy catheterization
- ✓ Softens and adjust to the anatomy
- ✓ Radiopaque for x-ray visibility

BUSTER Cat Catheter, side holes

- ✓ Soft and flexible catheter

Long term (0-7 days)

- Hospitalized patients and patients undergoing surgery
- The catheter has two lateral side holes and a rounded, atraumatically closed end
- Holes in the base plate allow suturing to the prepuce
- Luer lock on hub for easy connection to syringe or urine bag

BUSTER Barium Catheter

- ✓ Made of polyurethane and contains barium
- ✓ Radiopaque for x-ray visibility
- ✓ Polyurethane softens at body temperature increasing the biocompatibility. Hence, it is suitable for longer term placement and hospitalized patients