

# PERFUVET® TWIST

ROTATING INFUSION FIXATION SET



injection




surgery




orthopedics



hospitalization



capture



protection



hygiene



livestock



grooming

**génia**  
the veterinarian's preference.

## PERFUVET® TWIST

**NEW**


Anchored, rotating infusion fixation set that prevents the animal from tangling the tubing or blocking fluid flow. Ideal for fractious hospitalised patients.

- One base to attach to the cage equipped with double-sided tape for secure, fast attachment. Contains four holes to allow screwing into place.
- Two clasps that attach the infusion tubing to the base, allowing unlimited rotation.
- One 3 m length of sterile spiral PVC tubing with a Luer-Lock tip and a ring for attaching to the collar as needed.
- Perfuvet® Twist is easily attached to the cage; the base can be left on the cage and reused as often as necessary.
- The unlimited rotation prevents the animal from becoming tangled in the infusion set and blocking flow or requiring a time-consuming change of tubing.
- The spiral tubing allows ease of movement and prevents excessive pulling on the catheter.
- Saves time, improves safety of hospitalised animals, easy to use.
- Use on any size of animal (dogs, cats, rabbits...).

Sold as a kit (One fixed base, one clasp and three sets of single-use spiral tubing).


**GONE ARE THE DAYS OF ANIMALS GETTING TANGLED IN THEIR INFUSION SET, PERFUVET® TWIST, THE REVOLUTION IN PET HOSPITALISATION!**

Code	Model	Pack.
157450	Perfuvet® Twist kit (1 base, 2 clasps and 3 sets of single-use spiral tubing)	Unit
157451	Single-use spiral tubing	x 3



Watch the video version on our Youtube channel !

A TIME SAVER | EASY TO USE | SAFETY IN HOSPITALISATION



INJECTION