

SURGICRYL®

PGA

POLYGLYCOLIC ACID



Type	Braided and coated multifilament
Composition	Polyglycolic acid
Coating	Polycaprolactone and calcium stearate (1%)
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures



For quality
&
safety

smi

www.sutures.be

SURGICRYL[®]

PGA
POLYGLYCOLIC ACID

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

INDICATIONS General soft tissue closing and/or ligation
especially in General surgery
Skin closure
Gastrointestinal surgery
Gynaecology
Obstetrics
Plastic surgery
Urology
Ophthalmic surgery
Orthopaedics

BENEFITS

Excellent tissue support
Minimal tissue reaction
Predictable absorption rate
Excellent handling properties
High flexibility
High tensile strength
Excellent knot security
Soft passage through the tissues with minimal sawing, tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action
Guaranteed seal and product sterility

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery
Being absorbable not to use when an extended wound support is required.
Caution in patients with delayed wound healing (e.g. compromised immune system).

Without needle 

Metric	USP	Length		Order number
2	3/0	150 cm	<i>violet</i>	11200015
3	2/0	150 cm	<i>violet</i>	11300015
3.5	0	150 cm	<i>violet</i>	11350015
4	1	150 cm	<i>violet</i>	11400015
5	2	150 cm	<i>violet</i>	11500015
6	3&4	150 cm	<i>violet</i>	1160 0015
7	5	150 cm	<i>violet</i>	1170 0015
8	6	150 cm	<i>violet</i>	11800015
1.5	4/0	3 x 50 cm	<i>violet</i>	11150050
2	3/0	3 x 50 cm	<i>violet</i>	11200050
3	2/0	3 x 50 cm	<i>violet</i>	11300050
3.5	0	3 x 50 cm	<i>violet</i>	11350050
4	1	3 x 50 cm	<i>violet</i>	11400050
5	2	3 x 50 cm	<i>violet</i>	11500050
2	3/0	12 x 45 cm	<i>violet</i>	11200045
3	2/0	12 x 45 cm	<i>violet</i>	11300045
3.5	0	12 x 45 cm	<i>violet</i>	11350045
4	1	12 x 45 cm	<i>violet</i>	11400045
5	2	12 x 45 cm	<i>violet</i>	11500045
1.5	4/0	6 x 45 cm	<i>violet</i>	11150045-6
2	3/0	6 x 45 cm	<i>violet</i>	11200045-6
3	2/0	6 x 45 cm	<i>violet</i>	11300045-6
3.5	0	6 x 45 cm	<i>violet</i>	11350045-6
4	1	6 x 45 cm	<i>violet</i>	11400045-6
5	2	6 x 45 cm	<i>violet</i>	11500045-6
2	3/0	10 x 60 cm	<i>violet</i>	11200060
3	2/0	10 x 60 cm	<i>violet</i>	11300060
3.5	0	10 x 60 cm	<i>violet</i>	11350060
4	1	10 x 60 cm	<i>violet</i>	11400060
5	2	10 x 60 cm	<i>violet</i>	11500060
2	3/0	5 x 70 cm	<i>violet</i>	11200070
3	2/0	5 x 70 cm	<i>violet</i>	11300070
3.5	0	5 x 70 cm	<i>violet</i>	11350070
4	1	5 x 70 cm	<i>violet</i>	11400070
5	2	5 x 70 cm	<i>violet</i>	11500070







Without needle





Ligapack





Metric	USP	Length	Order number
1.5	4/0	250 cm	<i>violet</i> 11150025
2	3/0	250 cm	<i>violet</i> 11200025
3	2/0	250 cm	<i>violet</i> 11300025
3.5	0	250 cm	<i>violet</i> 11350025
4	1	250 cm	<i>violet</i> 11400025
5	2	250 cm	<i>violet</i> 11500025
6	3&4	250 cm	<i>violet</i> 11600025




Round bodied taper point




1/2 Circle






	Metric	USP	Length	Order number
HR 13 mm 	0.7	6/0	75 cm	<i>violet</i> 11070113
	1	5/0	75 cm	<i>violet</i> 11100113
	1.5	4/0	75 cm	<i>violet</i> 11150113
HR 17 mm 	0.7	6/0	75 cm	<i>violet</i> 11070117
	1	5/0	75 cm	<i>violet</i> 11100117
	1.5	4/0	75 cm	<i>violet</i> 11150117
	2	3/0	75 cm	<i>violet</i> 11200117
	3	2/0	75 cm	<i>violet</i> 11300117
HR 22 mm 	1.5	4/0	75 cm	<i>violet</i> 11150122
	2	3/0	75 cm	<i>violet</i> 11200122
	3	2/0	75 cm	<i>violet</i> 11300122
HR 26 mm 	1.5	4/0	75 cm	<i>violet</i> 11150126
	2	3/0	75 cm	<i>violet</i> 11200126
	2	3/0	75 cm	<i>undyed</i> 11200126U
	3	2/0	75 cm	<i>violet</i> 11300126
	3	2/0	75 cm	<i>undyed</i> 11300126U
	3.5	0	75 cm	<i>violet</i> 11350126
	4	1	75 cm	<i>violet</i> 11400126
HR 30 mm 	2	3/0	75 cm	<i>violet</i> 11200130
	3	2/0	75 cm	<i>violet</i> 11300130
	3	2/0	75 cm	<i>undyed</i> 11300130U
	3.5	0	75 cm	<i>violet</i> 11350130
	4	1	75 cm	<i>violet</i> 11400130
HR 36 mm 	2	3/0	75 cm	<i>violet</i> 11200136
	3	2/0	75 cm	<i>violet</i> 11300136
	3.5	0	90 cm	<i>violet</i> 11350136
	4	1	90 cm	<i>violet</i> 11400136




Round bodied taper point 		Metric	USP	Length	Order number
1/2 Circle 	HR 40 mm 	3	2/0	75 cm	<i>violet</i> 11300140
		3.5	0	90 cm	<i>violet</i> 11350140
		4	1	75 cm	<i>violet</i> 11400140-75
		4	1	90 cm	<i>violet</i> 11400140
		5	2	90 cm	<i>violet</i> 11500140
	HR 48 mm 	3	2/0	75 cm	<i>violet</i> 11300148
		3.5	0	90 cm	<i>violet</i> 11350148
		4	1	90 cm	<i>violet</i> 11400148
		5	2	90 cm	<i>violet</i> 11500148
		6	3&4	90 cm	<i>violet</i> 11600148








Round bodied taper point 		Metric	USP	Length	Order number
3/8 Circle 	DR 17 mm 	1.5	4/0	75 cm	<i>violet</i> 11150217
		2	3/0	75 cm	<i>violet</i> 11200217
	DR 25 mm 	3.5	0	75 cm	<i>violet</i> 11350225
		4	1	75 cm	<i>violet</i> 11400225




Round bodied taper point 		Metric	USP	Length	Order number
5/8 Circle 	FR 27 mm 	2	3/0	75 cm	<i>violet</i> 11200327
		3	2/0	75 cm	<i>violet</i> 11300327




Round bodied taper point 		Metric	USP	Length	Order number
J-needle 	JR 36 mm 	3	2/0	75 cm	<i>violet</i> 11303036
		3.5	0	75 cm	<i>violet</i> 11353036

Round bodied taper point 		Metric	USP	Length	Order number
5/8 Circle 	JR 36 mm  FR 27 mm 	3.5	0	75 cm	<i>violet</i> 11353036-0327
J-needle 					

Round bodied taper point 		Metric	USP	Length	Order number
1/2 Circle 	HR 50 mm 	4	1	150 cm	<i>violet</i> 11400950

Round bodied tapercut point 		Metric	USP	Length	Order number
1/2 Circle 	HRT 26 mm 	3	2/0	75 cm	<i>violet</i> 11300526
	HRT 30 mm 	3	2/0	75 cm	<i>violet</i> 11300530
		3.5	0	75 cm	<i>violet</i> 11350530
		4	1	75 cm	<i>violet</i> 11400530
	HRT 36 mm 	3	2/0	90 cm	<i>violet</i> 11300536
		3.5	0	90 cm	<i>violet</i> 11350536
		4	1	90 cm	<i>violet</i> 11400536
	HRT 40 mm 	3.5	0	90 cm	<i>violet</i> 11350540
		4	1	90 cm	<i>violet</i> 11400540
		5	2	90 cm	<i>violet</i> 11500540
	HRT 48 mm 	4	1	90 cm	<i>violet</i> 11400548





Round bodied tapercut point 		Metric	USP	Length	Order number
J-needle 	JRT 31 mm 	4	1	75 cm	<i>violet</i> 11403131









Reverse cutting 		Metric	USP	Length	Order number
3/8 Circle 	DS 12 mm 	0.7	6/0	45 cm	<i>violet</i> 11071512
		0.7	6/0	45 cm	<i>undyed</i> 11071512U
		1	5/0	45 cm	<i>violet</i> 11101512
		1.5	4/0	75 cm	<i>violet</i> 11151512




Reverse cutting






		Metric	USP	Length		Order number
3/8 Circle 	DS 12 mm Easy Pass	0.7	6/0	50 cm	<i>undyed</i>	11071512UP
		1	5/0	50 cm	<i>undyed</i>	11101512UP
		1.5	4/0	50 cm	<i>undyed</i>	11151512UP
	DS 16 mm	0.7	6/0	75 cm	<i>violet</i>	11071516
		1	5/0	75 cm	<i>violet</i>	11101516
		1	5/0	75 cm	<i>undyed</i>	11101516U
		1.5	4/0	75 cm	<i>violet</i>	11151516
		1.5	4/0	75 cm	<i>undyed</i>	11151516U
		2	3/0	75 cm	<i>violet</i>	11201516
	DS 16 mm Easy Pass	1.5	4/0	50 cm	<i>undyed</i>	11151516UP
	DS 19 mm	1	5/0	75 cm	<i>violet</i>	11101519
		1.5	4/0	75 cm	<i>violet</i>	11151519
		1.5	4/0	75 cm	<i>undyed</i>	11151519U
		2	3/0	75 cm	<i>violet</i>	11201519
		2	3/0	75 cm	<i>undyed</i>	11201519U
		3	2/0	75 cm	<i>violet</i>	11301519
	DS 19 mm Easy Pass	1.5	4/0	50 cm	<i>undyed</i>	11151519UP
		2	3/0	50 cm	<i>undyed</i>	11201519UP
		3	2/0	50 cm	<i>undyed</i>	11301519UP
	DS 24 mm	2	3/0	75 cm	<i>violet</i>	11201524
		2	3/0	75 cm	<i>undyed</i>	11201524U
		3	2/0	75 cm	<i>violet</i>	11301524
		3	2/0	75 cm	<i>undyed</i>	11301524U
		3.5	0	75 cm	<i>violet</i>	11351524
3.5		0	75 cm	<i>undyed</i>	11351524U	
DS 24 mm Easy Pass	2	3/0	50 cm	<i>undyed</i>	11201524UP	
	3	2/0	50 cm	<i>undyed</i>	11301524UP	
DS 30 mm	1.5	4/0	75 cm	<i>violet</i>	11151530	
	2	3/0	75 cm	<i>violet</i>	11201530	
	3	2/0	75 cm	<i>violet</i>	11301530	
	3	2/0	75 cm	<i>undyed</i>	11301530U	
	3.5	0	75 cm	<i>violet</i>	11351530	
	4	1	75 cm	<i>violet</i>	11401530	
	5	2	75 cm	<i>violet</i>	11501530	

Reverse cutting 		Metric	USP	Length	Color	Order number
3/8 Circle 	DS 36 mm 	3	2/0	75 cm	<i>violet</i>	11301536
		3.5	0	75 cm	<i>violet</i>	11351536
		4	1	75 cm	<i>violet</i>	11401536
	DS 40 mm 	5	2	75 cm	<i>violet</i>	11501536
		3	2/0	75 cm	<i>violet</i>	11301540
		3.5	0	75 cm	<i>violet</i>	11351540
		4	1	75 cm	<i>violet</i>	11401540
		5	2	75 cm	<i>violet</i>	11501540

Reverse cutting 		Metric	USP	Length	Color	Order number	
1/2 Circle 	HS 22 mm 	2	3/0	75 cm	<i>violet</i>	11201422	
		2	3/0	75 cm	<i>undyed</i>	11201422U	
		3	2/0	75 cm	<i>violet</i>	11301422	
		3	2/0	75 cm	<i>undyed</i>	11301422U	
	HS 26 mm 	1.5	4/0	75 cm	<i>violet</i>	11151426	
		2	3/0	75 cm	<i>violet</i>	11201426	
		3	2/0	75 cm	<i>violet</i>	11301426	
		3.5	0	75 cm	<i>violet</i>	11351426	
	HS 30 mm 	5	2	90 cm	<i>violet</i>	11501430	
		HS 36 mm 	3	2/0	90 cm	<i>violet</i>	11301436
			3.5	0	90 cm	<i>violet</i>	11351436
			4	1	90 cm	<i>violet</i>	11401436
	5		2	90 cm	<i>violet</i>	11501436	
	HS 40 mm 	3	2/0	90 cm	<i>violet</i>	11301440	
		3.5	0	90 cm	<i>violet</i>	11351440	
		4	1	90 cm	<i>violet</i>	11401440	
		5	2	90 cm	<i>violet</i>	11501440	
	HS 48 mm 	3.5	0	90 cm	<i>violet</i>	11351448	
4		1	90 cm	<i>violet</i>	11401448		
5		2	90 cm	<i>violet</i>	11501448		
6		3&4	90 cm	<i>violet</i>	11601448		

Cutting 		Metric	USP	Length		Order number
Straight 	GS 60 mm 	2	3/0	75 cm	<i>violet</i>	11202160
		2	3/0	75 cm	<i>undyed</i>	11202160U
		3	2/0	75 cm	<i>violet</i>	11302160

Cutting 		Metric	USP	Length		Order number
3/8 Circle 	NDS 16 mm 	1	5/0	45 cm	<i>violet</i>	11101716