



	Instant (0-10 minutes)	Flushing (0-60 minutes)	Short term (0-1 day)	Medium term (1-3 days)	Long term (0-7 days)
<b>Clinical Applications</b>	<ul style="list-style-type: none"> <li>Urine sample collection</li> <li>Check of passage in urethra</li> <li>Bladder tension relief</li> <li>Bladder lavage</li> </ul>	<ul style="list-style-type: none"> <li>Urine sample collection</li> <li>Check of passage in urethra</li> <li>Bladder tension relief</li> <li>Bladder lavage</li> <li>Flushing of urethra</li> </ul>	<ul style="list-style-type: none"> <li>Urine sample collection</li> <li>Check of passage in urethra</li> <li>Bladder tension relief</li> <li>Bladder lavage</li> </ul>	<ul style="list-style-type: none"> <li>Urine sample collection</li> <li>Check of passage in urethra</li> <li>Bladder tension relief</li> <li>Bladder lavage</li> <li>Catheterization for drainage purposes during surgery or hospitalization</li> </ul>	<ul style="list-style-type: none"> <li>Urine sample collection</li> <li>Bladder tension relief</li> <li>Bladder lavage</li> <li>Catheterization for drainage purposes during surgery or hospitalization</li> </ul>
<b>Recommended Catheters</b>	<p><b>BUSTER HC Urinary Catheter</b></p> <ul style="list-style-type: none"> <li>No need for gel</li> <li>Soft and comfortable, made of PVC</li> <li>Color-coded hub and label for easy size recognition</li> <li>Luer Lock on hub allows easy connection to syringe</li> </ul>	<p><b>BUSTER Dog Flushing Catheter</b></p> <ul style="list-style-type: none"> <li>More rigid catheter, made of polypropylene</li> <li>Luer Lock on hub allows easy connection to syringe or urine bag</li> <li>Atraumatic, rounded end hole for optimal flushing</li> </ul>	<p><b>BUSTER Dog Catheter</b></p> <ul style="list-style-type: none"> <li>More rigid catheter, made of polypropylene</li> <li>Luer Lock on hub allows easy connection to syringe or urine bag</li> <li>Atraumatic sides holes</li> </ul>	<p><b>BUSTER NC Urinary Catheter</b></p> <ul style="list-style-type: none"> <li>Soft and comfortable</li> <li>Colour-coded hub and label for easy size recognition</li> <li>Luer lock on hub allows easy connection to syringe or urine bag</li> </ul>	<p><b>BUSTER Foley with stylet and BUSTER Foley Plus without stylet</b></p> <ul style="list-style-type: none"> <li>Very soft silicone catheter with high biocompatibility</li> <li>Colour-coded hub and label for easy size recognition</li> <li>Urine drainage port fits KRUUSE urine bag</li> <li>Radiopaque for x-ray visibility</li> </ul>